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The screenshot displays the Philips software interface. On the left, there is a sidebar with patient information fields (Registration ID, Date of Birth, Gender) and a 'Plan' button. The main area shows a folder tree under 'Hospital'. A folder named '20XX' is selected, and a context menu is open over it. The menu items are: 'Unlock ExamCards Database', 'Copy', 'Select All', 'Import ExamCards...', 'Export ExamCards...', and 'Refresh'. The 'Unlock ExamCards Database' option is highlighted with a red box. At the bottom of the interface, a status bar displays the message: '15:56 ExamCards database has been locked and cannot be modified.'

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15:58 ExamCards database has been unlocked and can be modified.

Registration ID: 20XX\_101  
Date of Birth: 01/01/1999  
Gender: Male

EXAMINATION\_15..., 10/15/2024

Emotion\_training 00:46:29

- Survey
- EPI\_V1emotion EPI
- EPI\_V2emotion EPI
- EPI\_A1emotion EPI
- EPI\_A2emotion EPI
- MPRAGE 0.7iso sag

Accept

Summary	Geometry	Contrast	Motion	Dyn/Ang	Postproc	Offc/Ang	Coils	Conflicts
Angio / Contrast enh.	no				Total scan duration			10:10.0
Quantitative flow	no				Rel. SNR			1
Manual start	yes				Act. TR/TE (ms)			2000 / 30
Dynamic study	individual				Dyn. scan time			00:02.0
dyn scans				302	Time to k0			0.992
dyn scan times	shortest				ACQ voxel MPS (mm)			3.00 / 3.05 / 4.00
dummy scans	0				REC voxel MPS (mm)			3.00 / 3.00 / 4.00
immediate subtraction	no				Scan percentage (%)			98.4
fast next scan	no				Packages			1
synch. ext. device	yes				Min. slice gap (mm)			0
start at dyn.	1				EPI factor			63
interval (dyn)	1				Act. WFS (pix) / BW (Hz)			16.892 / 25.7
dyn stabilization	regular				BW in EPI freq. dir. (Hz)			2135.7
prospect. motion corr.	no				Min. WFS (pix) / Max. BW...			13.710 / 31.7
Keyhole	no				Min. TR/TE (ms)			1948 / 21
Arterial Spin labeling	no				Local torso SAR			< 60%
					Whole body SAR / level			< 0.8 W/kg / normal
					SED			< 0.5 kJ/kg
					Max B1+rms			1.26 uT
					PNS / level			74% / normal

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


EPI_V1emotion		Voxel	Tra
10:10		3.00 x 3.05 x 4.00	
Summary	<b>Geometry</b>	Contrast	Motion
IVIM SENSE: no			
k-t Acceleration: no			
Stacks: 1			
type: parallel			
slices: 32			
Slice thickness (mm): 4			
slice gap: user defined			
gap (mm): 0			
slice orientation: transverse			
fold-over direction: AP			
fat shift direction: P			
Minimum number of pa...: 1			
Slice scan order: FH			
Large table movement: no			
PlanAlign: no			
REST slabs: 0			
Patient position: head first			
Patient body position: head first			
Patient orientation: supine			
Patient body orientation: supine			

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